



labocool

Laboratory Fridges

**Cold storage protection
for your valuable
samples & materials**



**Single &
Double
Door**

- LED Illumination
- Automatic Defrost
- Magnetic Door Gasket
- Adjustable Shelves
- ISO9001
- ISO14001
- 2 Year Warranty
- Double Tempered Glass
- Best in Class Components
- Manual Thermostat Preset at 3.5°C
- Soft Closing Doors
- Plastic Coated Shelves



LABOTEC

Your Partner in Science & Technology

Head Office (011) 315 5434 • Cape Town (021) 531 7660 • Durban (031) 566 4870
sales@labotec.co.za • www.labotec.co.za

The Labocool fridges from Labotec have an operating range from +1°C to +8°C. They are manufactured in accordance with ISO 9001 (Quality management systems certified) and ISO 14001 (Environmental Management system compliant) and come with a two year warranty for added peace of mind.

The range is ideally suited for general storage and typically can be used for reagents and kits, whole blood chromatography apparatus, test samples, biological and culture media and more. The glass door is double tempered for safety, and the magnetic door gasket promotes stable temperature within the cabinet.

Space-saver Fridge 224 Lt



Dimensions		
(mm)	External	Internal
width	586	495
depth	605	500
height	1340	905
weight	83kg Net	
Foot Print Area	0.35m ²	
Gross Volume	224Lt	
Number of shelves (STD)	3	
Shelf Size	w472 x d458(mm): 0.2m ²	

Space-saver Fridge 627 Lt



Dimensions		
(mm)	External	Internal
width	710	625
depth	730*	644
height	2040	1508
weight	133kg Net	
Foot Print Area	0.61m ²	
Gross Volume	627Lt	
Number of shelves (STD)	5	
Shelf Size	w595 x d542(mm): 0.32m ²	

*Excluding 40mm handle & 100mm spacer

Double Door Fridge FV 1000



Dimensions		
(mm)	External	Internal
width	1025	941
depth	731*	644
height	2045	1508 1396
weight	133kg Net	
*Excluding 100mm spacer		
Foot Print Area		0.75m ²
Gross Volume		936Lt
Number of shelves (STD)		5x5
Shelf Size	w458 x d543(mm): 0.25m ²	

Shelving is adjustable so it can be moved around for optimal usage of the cabinets' capacity.

The shelves are also plastic coated to ensure that they stay hygienic and rust-free. An automatic defrost minimises any trace amounts of frost build-up.

Optional Extras



Electronic Voltage Protection



Lock



Digital Thermometer

- **Magnetic Door Gasket**
Promotes Stable temperature within cabinet
- **Adjustable Shelves**
Ensures optimum usage of cabinet capacity
- **Plastic Coated Shelves**
Rust free, hygienic, long lasting
- **2 Year Warranty**
Long term peace of mind
- **Double Tempered Glass**
Adhering to SA Safety requirements
- **Best in Class Components**
Peace of mind with a quality product
- **Automatic Defrost**
Minimizes trace amounts of frost build up
- **ISO9001**
Quality Management Systems Certified
- **ISO14001**

Environmental Management Systems Compliant
(Entire Manufacturing process is performed in an environmentally friendly manner)



labocool